

THE #1 PRESCRIBED CUSTOM ANTERIOR BITE SPLINT IN THE WORLD!



BiteSoft™ Anterior Splint

**Easily Manage Bruxism
Protect Your Patients Restorative Investment
Increase Patient Compliance
Save Valuable Chairtime**

COVERS ONLY 6 ANTERIOR TEETH

Provides Ideal Retention and is Extremely Comfortable to Wear
Worn Only at Night, Provides Minimal Risk for Occlusal Change

LINGUAL RAMP

Ensures Mandible Stability and Relaxation of Tense Jaw Muscles to Help Reduce Headaches and Jaw Pain



TWO PREMIUM MATERIALS AVAILABLE:
Dual Laminate and Thermoplastic Lining

UNIQUE PRODUCT Authenticity Code (PAC)

Indicates the Legitimacy of Each BiteSoft Splint

CENTRIC OCCLUSION TABLE

Contact of Only 4 Lower Incisor Teeth in Centric Occlusion Reduces Forces of Bruxism by up to 70%



BITESOFT EDUCATION VIDEO
SCAN with your Smart Phone to view or online at www.tridentlab.com/BiteSoftVideo



Request Your FREE Patient Education Demo Model!



CALL TODAY FOR A CASE PICK-UP!
(800) 221-4831

www.tridentlab.com



12000 Aviation Boulevard Hawthorne, CA 90250 • Fax (310) 915-7171

PROVIDING INNOVATIVE DENTAL PRODUCTS AND COMMITMENT TO EXCELLENCE SINCE 1988